



**EST.1991**



**LEODIS**  
**STAIRLIFTS**



**HELPING YOU WITH THE UPS AND DOWNS IN LIFE**



**STRAIGHT STAIRLIFTS**

# Are you having problems or struggling to get up and down your stairs?

Loss of mobility throughout the home can be a frustrating and disheartening experience. The decision of whether to stay in your home or move house is a very difficult choice. Having a stairlift fitted could be the ultimate solution to a difficult problem allowing you to stay with the safety and security of your own home.

Established since 1991, **LEODIS** STAIRLIFTS have helped people regain their independence in their own home, offering them choice from leading manufacturers, impartial advice and a personal service, as expected of a family run company, that is second to none.


At **LEODIS** STAIRLIFTS our aim is to guide you through the options available to you to suit your specific requirements and make your experience with us a smooth and easy process. Our service and support is designed to provide you with confidence and complete peace of mind.

## Straight Stairlifts

Quality, safety and reliability are foremost when buying a stairlift and here at **LEODIS** STAIRLIFTS we only supply the best. As agents for the leading manufacturers, we can offer you impartial advice on the straight stairlift that is right for your needs both now and in the future.

Buying any stairlift can be confusing so in this brochure we aim to explain the options available to you, enabling you to make a more informed choice. A free, no obligation, home survey is an essential part of the **LEODIS** service, here we will advise you which of our products best suits your needs and your budget. We will measure your staircase and discuss if any alterations need to be made i.e. moving a radiator and, in areas where access may be obstructed, if you need a hinged/sliding track.

**LEODIS** STAIRLIFTS supply and install both new and reconditioned straight stairlifts. Our stairlifts are available for installation within days of your order with us and we can, subject to stock, offer NEXT DAY installations.



# Guiding you through the process of buying a stairlift from **LEODIS STAIRLIFTS**

Our service and support at **LEODIS STAIRLIFTS** is designed to provide you with confidence and complete peace of mind. We want to make your experience with us as hassle free as possible.



1. During our initial conversation we will arrange your **FREE**, no obligation home survey.
2. At your **free** home survey we will:
  - Discuss the options available to you, as agents for the leading manufacturers, we can offer you **impartial** advice.
  - Show you photographs of our stairlifts in situ.
  - Measure / survey your staircase.
  - leave you with a written quotation and our selling pack.

We will then leave you alone.

We pride ourselves on our **NO hassle / NO call** back service.

3. Should you wish to go ahead, we will arrange a mutually convenient date and time for your stairlift installation. This can be often the **NEXT DAY**.

You will not be asked for a deposit.

4. Your installation will be carried out by one of our experienced, friendly stairlift engineers. He will install your stairlift and demonstrate how your stairlift will work, answering any questions/queries you may have. Your home will be left in a clean and tidy fashion.

Only when you are entirely satisfied do you pay.

Should you need to contact **LEODIS STAIRLIFTS** you can rest assured your call will be answered in person by an engineer  
24 hours a day, 7 days a week  
(including evenings, weekends and bank holidays)

# Handicare / Minivator simplicity 950 straight stairlift

The ultra slim, simplicity 950 straight stairlift, supplied and fitted by **LEODIS** STAIRLIFTS, is the simple solution for those who want a dependable stairlift but do not require any optional extras such as a powered swivel and / or footplate. Fitted with an easy to swivel seat for your safe exit and entry onto the stairlift, the simplicity 950 seat comes in a neutral (sand) shade and has a clip seatbelt as standard. Available new or reconditioned for NEXT day delivery.

The beauty of the simplicity range is that the batteries will charge at any point on the track so the stairlift can be positioned anywhere when not in use. The simplicity 950 comes with 2 remotes as standard so can be called from any position on the track.

The simplicity 950 can have a sliding or manual hinged track if required.

Please see details on simplicity 950+ page.



## How do I operate my stairlift?

The unique toggle is standard on all Handicare / Minivator stairlifts. It can be situated on the end of either arm rest and can be easily operated by hand or arm.



Call / Send remotes are supplied as standard with our simplicity range of stairlifts, allowing you to call and send your stairlift without being seated on it.



### Useful information

Footplate to seat: 390-560mm  
Width of seat at widest: 480mm  
Min. folded intrusion: 285mm  
Max. user weight: 140kg/22st  
Max. user weight with hinged / sliding track: 120kg / 19st

## Handicare / Minivator simplicity 950+ straight stairlift

The Handicare / Minivator 950+ straight stairlift offers the same reliability and elegant, clean lines of the 950 but comes with a choice of 6 seat colour options to suit any decor and the option for a powered swivel seat and/or powered footplate. The inertia seat belt fastens either to the side or centrally for your comfort. This stairlift is available new and reconditioned for next day delivery from **LEODIS** STAIRLIFTS.

The Simplicity 950+ comes with 2 remotes as standard so can be called from its standing position on the track. The ultra slim design fits most straight staircases whilst still allowing plenty of access to the stairs for others. As with the Simplicity 950, for your convenience, the stairlift will charge at any point on the track.




### Do you have restricted space in your home?

For areas with restricted access a hinged track may be appropriate. The track is moved manually and a gas assisted hinge makes operation easy, however a powered slide track is also an option.

---

Handicare / Minivator have pioneered a solution that removes the need for a hinge. The sliding track automatically travels up the staircase at the same time as your stairlift, leaving the doorway / hall obstruction free. Ask your **Leodis** engineer for more information.

---



#### Useful information


Footplate to seat: 430-540mm

Width of seat at widest: 460mm

Min. folded intrusion: 370mm

Max. user weight: 140kg/22st (120kg/19st with hinged/sliding track)

#### seat colour options



# Handicare / Minivator Perch / Stand Straight Stairlift

It is our aim at **LEODIS** STAIRLIFTS to find you the stairlift that suits your needs, for the majority of people that involves a 'sit on' stairlift. If you have painful knee or hip joints or restricted movement there may be another option for you.



The Minivator / Handicare Perch seat allows you to remain in a standing position whilst accessing your stairs. The seat is available with the Simplicity 950+ and comes with the standard features of that model including an inertia seat belt for your safety.

The seat comes in a choice of 6 colours to suit any décor.



#### Useful information

Footplate to seat: max 870mm

Width of seat at widest: 460mm

Max. user weight: 140kg/22st

For further measurements

please ask you **LEODIS** STAIRLIFT engineer.

---

## Handicare / Minivator 1000 outdoor straight stairlift

At **LEODIS** STAIRLIFTS we understand the problems associated with the loss of mobility around your home. For many people steps out into the garden can also be problematic, making them feel they are trapped within their own home.

**LEODIS** STAIRLIFTS may be able to offer you the perfect solution.

The Minivator 1000 outdoor stairlift is equipped with the standard features expected of a Handicare / Minivator stairlift however is robust enough for outdoor use. The upholstered seat has a waterproof cover to protect it from the elements.



#### Useful information

Footplate to seat: 450-570mm

Width of seat at widest: 440mm

(+50mm with arm spacers)

Min. folded intrusion: 387mm

Max. user weight: 160kg/25st

To prevent unauthorised use of the lift simply isolate the key switch until you wish to use your stairlift again.



## Brooks 120/130 slim-line straight stairlift

Brooks are among the leading manufacturers in stairlifts. The Brooks 120/130 Slim-line straight stairlift offers you the reliability and quality as expected of a leading manufacturer. Whilst looking at home on any straight staircase the Brooks 120/130 Slim-line still offers plenty of access to the stairs when folded. An easy-to-use paddle control operates this smooth, comfortable lift whilst the inertia seat belt is fitted for your safety.

After taking you to your destination the lift is designed to rest on its charge points, situated either end of the track, so is then ready for your next use, it comes with 2 remotes as standard. The Brooks 120/130 is available for NEXT DAY installation from **LEODIS** STAIRLIFTS.



### A solution for small spaces

If your track goes across a doorway or an access area then you may require a hinged track. The track is easily operated manually, however a powered option is also available.

(Manual and powered hinged tracks are optional extras).

Ask your **LEODIS** engineer your for further details.



### How do I operate my stairlift?



Simply push the paddle switch in the direction in which you wish to travel and your stairlift will take you quietly and effortlessly to where you want to go.



## Service and Support

At **LEODIS STAIRLIFTS** we aim to provide a first class service and support that is second to none, giving you and your family complete peace of mind.

We pride ourselves on the quality of our products and on all of our installations and whilst our stairlifts are completely reliable, from time to time, you may need our assistance or to speak to one of our engineers. By being a **LEODIS STAIRLIFTS** customer you are guaranteed:

- A locally based company that is contactable 24 hours of every day – important should you or your family have any questions or need a call out.
- Every call you make to us is answered in person by an engineer – YOU are our priority.
- 24 hour call out cover for all OUR customers, 7 days a week (including evenings, weekends & Bank Holidays).
- 2 FREE service visits on new models / 1 FREE service visit on reconditioned models.

## Warranties

**LEODIS** new straight stairlifts come with up to 10 years warranty and 2 FREE service visits. Our reconditioned straight stairlifts come with 12 months parts and labour warranty and include 1 FREE service visit.

## Safety as Standard

Your safety is of paramount importance to **LEODIS STAIRLIFTS**. The straight stairlifts that we supply and install have standard safety features, as expected of the leading manufacturers.

A commonly asked question is what happens if there is a power cut? All our stairlifts are battery operated and will still work for a limited time if you have a power cut.

All of our stairlifts are fitted with seatbelts as standard for your safety. For your comfort and ease of use they are easy to release.



Should your stairlift encounter an obstruction on the staircase the sensors on the footplate will immediately halt your stairlift. To release, operate the stairlift in the opposite direction and remove the object before continuing on your journey.



All **LEODIS STAIRLIFTS** are fitted with a lockable isolation key to prevent unauthorised use of your stairlift.



## Your Questions Answered

Hopefully by reading through our brochure we have answered any questions you may have, however we have answered some frequently asked questions below.

### How long will I have to wait for my stairlift?

Once you have placed an order with **LEODIS STAIRLIFTS** your stairlift can often be installed within days and, in some cases, as quickly as the very **NEXT** day.

### How long will it take to install my straight stairlift?

Our experienced, friendly engineer will fit your straight stairlift in approximately 2-3 hours.

### Will my stairlift be attached to the wall?

Many people think that stairlifts are fixed to the wall however this is not the case. The stairlift will be attached directly to the staircase itself. During your free home survey our **LEODIS** engineer will advise you on the most convenient side of the stairs for your stairlift and if anything ie radiator may need to be moved.

### What happens if there is a power cut?

Your stairlift will still work for a limited time if you have a power cut. Your stairlift is battery operated and these re-charge when your stairlift is left on a charge point top or bottom of the track or, in the case of the Simplicity 950, anywhere on the track. Your stairlift will alert you by bleeping when it is not left on a charge point. If your stairlift batteries go flat they may need to be replaced.

### Do you only sell new straight stairlifts?

At **LEODIS STAIRLIFTS** we can offer both new and reconditioned straight stairlifts, both available for NEXT day delivery in most cases. As agents for the leading manufacturers we can advise and install the product that best suits YOUR requirements.

### What do I do when I need a new carpet?

The best time to replace your carpet is before installation however if this is not possible **LEODIS STAIRLIFTS** can remove then reinstall your stairlift, a cost will be incurred for this service.

If you have any further questions or queries **LEODIS STAIRLIFTS** will be happy to answer your question in person. Call us NOW on **0800 542 6045**.



EST.1991

# LEODIS STAIRLIFTS

At **LEODIS STAIRLIFTS** we know that finding the solution for you and your home is the most crucial part of the service so we use our wealth of knowledge and experience to find the product that best suits YOUR specific needs.

Call us now on

**0800 542 6045**

## CONTACT US

Email: [sales@leodisstairlifts.co.uk](mailto:sales@leodisstairlifts.co.uk)  
[WWW.leodisstairlifts.co.uk](http://WWW.leodisstairlifts.co.uk)

